

CARE

RINSE THE PUMP IN WARM WATER AND TREAT THE SOFT SLEEVE AND CYLINDER WITH A SPECIAL TOY CLEANER BEFORE AND AFTER EACH APPLICATION.

STORAGE

DIRECT SUNLIGHT IMPACT CAN DAMAGE THE PUMP. KEEP THE DEVICE IN A DARK PLACE, OUT OF THE REACH OF CHILDREN, IN THE ORIGINAL PACKAGING.

USAGE

BEFORE STARTING TRAINING, IT IS RECOMMENDED TO CONDUCT A TRIAL SESSION LASTING 5 MINUTES. THIS WILL ALLOW YOU TO UNDERSTAND HOW THE DEVICE WORKS, AND HOW IT SHOULD BE HANDLED.

IT IS RECOMMENDED TO APPLY A LUBRICANT TO THE GENITAL ORGAN BEFORE USAGE. USING A PUMP WITH A LUBRICANT WILL HELP TO AVOID POSSIBLE DAMAGE TO THE PENIS. ALSO, IT WILL MOISTURIZE THE SKIN OF THE ORGAN AND RELIEVE IT FROM DRYNESS. THE LACK OF LUBRICANT OR USE OF UNDEDICATED STUFF (CREAM, VASELIN, OIL, FAT-CONTAINING MATERIALS) DESTROY THE MATERIAL OF THE SLEEVE.

THE MAIN FUNCTIONS OF THE VACUUM PUMP:

1. TRAINING TO INCREASE THE PENIS AND ERECTILE STRENGTH.
2. MAINTAINING AN ERECTION IN MEN WITH ERECTILE DYSFUNCTION.
3. PREVENTION OF PREMATURE EJACULATION WHEN USING A PUMP BEFORE SEX.
4. OBTAINING SEXUAL SATISFACTION AND PREVENTION OF DECREASED LIBIDO IN THE ABSENCE OF A SEXUAL PARTNER FOR A LONG TIME.

MANUAL

DISCOVERY

THE SCHEME OF USING A VACUUM PUMP FOR PENIS ENLARGEMENT

1. BEFORE STARTING TO PUMP AIR, IT IS NECESSARY TO PUT AN ERECTILE RING ON THE CYLINDER.
 2. START CREATING A VACUUM BY PUMPING AIR OUT OF THE CYLINDER.
 3. WHEN THE ORGAN REACHES AN ERECTILE STATE – MOVE THE RING TO THE BASE OF THE PENIS. THIS WILL HELP TO FIX THE RESULT, BLOCKING THE OUTFLOW OF BLOOD.*
 4. MAINTAIN THE DURATION OF THE TRAINING, THEN PRESS THE PRESSURE RELIEF VALVE.
 5. REMOVE THE ERECT AND ENLARGED PENIS FROM THE CYLINDER.
- * - A SILICONE RING TO MAINTAIN AN ERECTION IS NECESSARY IF THE USER HAS ERECTILE DYSFUNCTION.

MANUAL

DISCOVERY

ADDITIONAL INFORMATION

THE PUMP SHOULD FIT SNUGLY TO THE INGUINAL AREA ABOVE THE SCROTUM. THE USAGE OF THE PUMP SHOULD NOT CAUSE DISCOMFORT OR PAINFUL SENSATIONS. IF YOU EXPERIENCE PAIN OR DISCOMFORT, DRAIN THE AIR FROM THE PUMP OR STOP USING IT.

FOR BEST RESULTS, IT IS RECOMMENDED TO USE THE PUMP FROM 3 TO 5 TIMES A WEEK. EACH TRAINING SESSION SHOULD CONSIST OF SEVERAL APPROACHES.

THE AVERAGE DURATION OF PUMPING CAN BE UP TO 5 MINUTES WHILE MAINTAINING A COMFORTABLE FEELING. TO ACHIEVE THE BEST RESULT, IT IS NECESSARY TO ALTERNATE THE PUMPING OF AIR.

IT IS NOT RECOMMENDED TO EXCEED THE TIME OF USING THE PUMP FOR MORE THAN 20 MINUTES.

THERE ARE A NUMBER OF CONTRAINDICATIONS. CONSULT YOUR DOCTOR BEFORE USE.